Dental Implant and Abutment Instructions For Use

Indications Dental Implants

Intra-Lock® Implants have been designed to restore partially or fully edentulous patients. The implants have been designed to be used in either the mandible or the maxilla and to support removable or fixed prostheses, from single tooth replacement to full arch reconstruction. Intra-Lock® Implants are intended for immediate function on single tooth and/or multiple tooth applications when good primary stability is achieved, with appropriate occlusal loading, in order to restore normal teeth functions.

Indications MILO™ Dental Implant Systems

MILO™ Dental Implants are indicated for long-term maxillary and mandibular tissue-supported denture stabilization. They are also indicated for the rehabilitation of single or maxillary lateral incisors and mandibular lateral and central incisors. Multiple implants may be restored after a period of time or placed in immediate function.

Indications Mini Drive-Lock™ Dental Implant System

Mini Drive-Lock™ Dental Implants are intended for use as a self-tapping titanium screw for transitional or intra-bony long-term applications.

Mini Drive-Lock™ Dental Implants are also indicated for long-term maxillary and mandibular tissue-supported denture stabilization. Multiple implants should be used and may be restored after a period of time or placed in immediate function.

Indications Dental Implants Blossom

Intra-Lock® Implants have been designed to restore partially or fully edentulous patients. The implants have been designed to be used in either the mandible or the maxilla and to support removable or fixed prostheses, from single tooth replacement to full arch reconstruction. Intra-Lock® Implants are intended for immediate function on single tooth and/or multiple tooth applications when good primary stability is achieved, with appropriate occlusal loading, in order to restore normal teeth functions.

Indications OP Dental Implants

For 3.0 mm - Intra-Lock® OP Dental Implants are indicated for long-term maxillary and mandibular tissue supported denture stabilization. They are also indicated for rehabilitation of single or maxillary lateral incisors and mandibular lateral and central incisors. Multiple implants may be restored after a period of delayed loading or placed in immediate function when good primary stability is achieved with appropriate occlusal loading in order to restore normal teeth function.

For 3.75mm, 4.0mm & 4.75mm Intra-Lock OP Dental Implants

Intra-Lock® OP Dental Implants have been designed to restore partially or fully edentulous patients. The implants have been designed to be used in either the mandible or the maxilla and to support removable or fixed prostheses, from single tooth replacement to full arch reconstruction. Intra-Lock® Implants are intended for immediate function on single tooth

and/or multiple tooth applications when good primary stability is achieved, with appropriate occlusal loading, in order to restore normal teeth functions.

Indications Dental Abutments

Intra-Lock® abutments are intended for use with Intra-Lock® dental implants to support a prosthetic device in partially or fully edentulous patients. The abutments may be used in single and/or multiple tooth application in the mandible or maxilla.

Contraindications

Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies or patients who have had myocardial infarction within the last 12 months.

Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or non-compliance to oral hygiene procedures, or patients who cannot maintain oral hygiene procedures if implants are placed.

Tobacco usage increases the occurrence of complications and failures.

Material

Intra-Lock® implants and abutments are manufactured from Titanium 6AL-4V ELI Alloy, ASTM F136 or Commercially Pure Titanium, ASTM F-67. Intra-Lock® manufactures a temporary abutment from PEEK and Titanium 6AL-4V ELI.

How Supplied

Intra-Lock® dental implants are provided sterile (by gamma radiation) and are intended for single use only. Packaged implants are suspended on a titanium or polycarbonate support ring within a clear vial. This vial is placed into a Seal Pac™ plastic vial with a tamper evident seal, which provides an additional environmental barrier.

The label on the package provides the lot number, product name, catalog reference number and expiration date. To ensure sterility, dental implants must be used before the end of the expiration date indicated on the outer package label.

Never reuse, reclean or resterilize a dental implant. These activities can adversely affect implant materials and alter the surface characteristics, which may result in poor function and implant failure.

Prior to using Intra-Lock® dental implants, inspect the package and labeling for integrity. If the device is opened, damaged or contaminated in any way, DO NOT USE.

Intra-Lock® dental abutments are provided non-sterile and intended for single use only. Abutments should be sterilized prior to use by the instructions provided.

Prior to using Intra-Lock® dental abutments, inspect the package and labeling for integrity. If the device is opened, damaged or contaminated in any way, it must not be used.

Sterilization of Abutments

STERILE 

Abutments may be sterilized using a full cycle pre-vacuum steam sterilization at a temperature of 132°C for an exposure time of 3 minutes and 5 minute drying time.

Preoperative Treatment Planning

Proper patient selection is important. A comprehensive patient interview and medical/dental history must be taken. A complete oral examination should then be conducted. Head and neck examination is followed by a thorough oral examination. The use of magnification is strongly encouraged as an adjunct to all procedures.

Oral inspection includes palpation and the proper radiographic protocol(s). This may include periapicals, panorex and tomograms. Palpation of the ridges is also required and the use of intra-oral probes for tissue thickness is recommended.

The diagnostic procedures will give the dentist an appreciation for the tissue quality and thickness, ridge morphology for the type and size of the implants that might be required. Measurements for implant size can be estimated utilizing radiographs, templates, calipers and millimeter rulers.

Treatment planning should also take into consideration prosthetic biomechanics, occlusion and occlusal load. Fracture due to excessive load or metal fatigue can occur on the implant body or its prosthetic components if this aspect of planning is inadequate.

In overdenture cases four or more implants should be utilized for maxillary or mandibular tissue-supported denture stabilization.

When fixed prosthetics are utilized in single stage surgical procedures, implants may be loaded immediately following insertion provided at least four implants are placed and are splinted with a bar. These implants should be placed principally in the anterior mandible, between the mental foramina, where good initial stability of the implants can most often be achieved.

Surgical Asepsis

As with all surgical procedures, the operatory field should be maintained with sterile coverings (light handles, chair controls, bracket tray, and all instruments and components). Barrier technology, sterile solutions and sprays, sterile coverings, and proper autoclaving techniques must be employed as indicated.

Important Points to Remember

- The Conic Implants have a tapered thread body and are designed to fit the drilled depth.
- The final seating of the implant(s) should be achieved by the incremental turns of the Ratchet Wrench.
- Radiograph(s) should confirm proper depth, seating, orientation and placement of the implant(s).
- If a removable prosthesis is used near the implant site, in process of healing, it should be generously relieved and a soft tissue conditioner reline material placed.

Postoperative Care

Cold packs are recommended for the first 24 hours. Analgesics/Antibiotics may be prescribed at the discretion of the practitioner. The patient should be advised to favor the

opposite side of the mouth, maintain a soft diet and avoid hot liquids.

Warning

The Intra-Lock® Dental Implant System has not been evaluated for safety and compatibility in the MR environment. The Intra-Lock® Dental Implant System has not been tested for heating or migration in the MR environment.

Warnings

- Dental implant surgery is a complex dental procedure. Appropriate and adequate training in all phases of implant procedures and proper technique is strongly recommended prior to implant use.
- Improper patient selection, diagnosis, treatment planning or technique can result in implant failure and/or loss of supportive bone.
- Care must be taken if performing electrosurgery around a dental implant. Electrosurgery generates heat, which can be conducted through a metallic implant and cause damage to surrounding tissue or bone.
- The use of small diameter implants and angled abutments in the posterior region of the mouth is not recommended due to possible failure of the implant.
- The external surface of Intra-Lock® Dental Implants should only come in contact with titanium surfaced instruments.
- Prior to the placement of all dental implants and dental abutments, appropriate and adequate training in all phases of implant and abutment procedures is strongly recommended.
- A comprehensive patient interview, medical/dental history and complete oral examination should be conducted.
- Diagnostic radiographs and mounted study models, if appropriate should also be obtained.


Warning

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist or physician.

Intra-Lock® International Inc.
6560 West Rogers Circle, Bldg. 24
Boca Raton FL 33487 – USA

www.intra-lock.com

Intra-Lock System Europa, S.p.A.
I-84100 Salerno
+ 39 089 233 045

IFU-DIS (EN) (07/2018)

Dental Drill Instructions For Use

Indications

Intra-Lock® dental drills are for use to perforate the tissue or bone and create the osteotomy, or hole, where Intra-Lock® dental implant(s) will be placed.

Material

Dental drills are manufactured from Stainless Steel.

How Supplied

Sterile Dental Drills

Dental drills that are provided sterile (by gamma radiation) and are intended for single patient use only.

The label on the package provides the lot number, product description, catalog reference number and expiration date. To ensure sterility, dental drills must be used before the end of the expiration date indicated on the outer package label.

Prior to using dental drills, inspect the package and labeling for integrity. If the package has been opened or damaged or the device has been contaminated or exceeded the expiration date before the first use, return dental drill for replacement.

Non-Sterile Dental Drills

Non-sterile drills MUST be steam sterilized prior to use. See Steam sterilization instructions. The label on the package provides the lot number, product description and catalog reference number.

Prior to using dental drills, inspect the package and labeling for integrity. If the new package has been opened or damaged follow the instructions for steam sterilization.

Steam Sterilization

Dental drills may be sterilized using a full cycle pre-vacuum steam sterilization at a temperature of 132°C (270°F) for an exposure time of 4 minutes and a dry time of 30 minutes.

Warnings

- Steam Sterilization instructions have not been validated to reprocess dental drills for use with multiple patients. After use on any patient, the dental drill may include extensive tissue and/or bone debris and therefore require more extensive cleaning and sterilization techniques.
- If dental drills are reused, the user or reprocessing facility must validate the parameters they utilize, which must include an inspection for potential wear and damage incurred to the dental drill from the previous use.

Drilling Technique

Warnings

- Osseous drilling can generate detrimental heat that may cause thermal bone necrosis and subsequent failure of the dental implant to osseointegrate. It is critically important that all bone-cutting procedures are CONDUCTED AT SLOW SPEEDS. Profuse, internal and/or external irrigation is mandatory.
- The slow-speed, highly irrigated drilling procedure is conducted while angling the dental drill such that the direction of the drill bisects the ridge. The dental drill should also be held vertically, avoiding a mesial or distal cant.
- Depth gauge/alignment components can be periodically inserted into the osteotomy site to monitor the angle of penetration. Successively larger drills are used until the desired diameter is achieved.
- Drilling is performed with a precise, up and down pumping action. The drill angle is maintained in order to preserve the concentricity of the hole, while the pumping action allows for incremental depth penetration and periodic cleansing of the flutes.
- It is recommended to drill a Pilot Hole when the bone is of high density using a "Lancer" Pilot drill extending to a depth of 8 mm.


Warning

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist or physician.

Intra-Lock® International Inc.
6560 West Rogers Circle, Bldg. 24
Boca Raton FL 33487 – USA

www.intra-lock.com

Intra-Lock System Europa, S.p.A.
I-84100 Salerno
+ 39 089 233 045

0086

IFU-DR (EN) (07/2018)

ReOss® Synthetic Biomaterial Instructions for Use

Indications

ReOss® Powder and Putty is indicated for filling and/or augmenting intraoral/maxillofacial osseous defects, such as intrabony periodontal osseous defects, furcation defects, augmentation of bony defects of the alveolar ridge, filling of tooth extraction sites and sinus elevation grafting.

Contraindications

Standard protocols are to be followed when taking a comprehensive medical/dental history.

Uncontrolled diabetes, severe vascular or neurological disease, blood dyscrasias, existing acute or chronic infections, especially at the site of use, metabolic or systemic bone or calcium disorders, severe degenerative bone disease.

History of radiation to the treatment area, history of drug/alcohol abuse, history of allergic reactions to PLGA products

Renal compromised patients, uncooperative patient who will not or cannot follow postoperative instructions, including individuals who abuse drugs and/or alcohol, pediatric patients.

How Supplied

ReOss is intended for single use only. Packaged in a glass bottle with a bromide stopper and a tamper evident aluminum crimped seal.

The label on the package provides the lot number, product name, catalog reference number and expiration date. To ensure sterility, ReOss must be used before the end of the expiration date indicated on the outer package label.

Never reuse or resterilize ReOss. This activity can adversely affect materials, which may result in poor function and failure.

Prior to using ReOss, inspect the package and labeling for integrity. If the device is opened, damaged or contaminated in any way, DO NOT USE.

Material

ReOss® is a hydrophilic, highly porous, resorbable, synthetic copolymer permeated with osteoconductive particles of Hydroxyapatite. It is configured as a multi-pore three - dimensional scaffold that creates an environment suitable for the infusion of blood and osteoblasts and is engineered to integrate with the physiochemical state of bone tissue.

Pre-Processing

ReOss® guidelines are to be used in conjunction with established surgical techniques. Standard procedures for treatment of bone defects involving bone grafting and internal fixation should be followed. The results of bone grafting procedures are highly variable. When selecting the bone grafting material and the surgical technique to be utilized, the following factors should be taken into consideration:

- Bone quality and location of the defect;
- Anticipated load factor;
- Age of the patient;

- Suitable blood supply at graft site;
- Proper graft apposition;
- Complete coverage of the graft material to maintain material in situ.

For best results, care should be exercised to ensure the correct graft material is selected for the intended application.

Radiographic imaging of the defect is recommended to assess the location and extent of the defect to determine the optimum type and amount of ReOss® to be used.

Standard aseptic technique must be maintained to minimize the risk of post-operative complications.

Standard postoperative regimen as used when soft tissue and bone augmentation procedures should be followed. All cases must be treated in a manner consistent with individual patient profile, nature and extent of procedure.

Precautions

ReOss® is not intended for load-bearing applications.

The effects of ReOss® on pediatric patients is unknown.

⚠️ Warnings

- Pre-treatment, treatment and post-treatment protocols should be followed as per standard treatment procedures.
- Follow established protocols, including pre-and post-medication regimens prescribed to avoid infection, when/where necessary.
- Avoid overfilling and do not disturb the material until hardened.
- Soft tissue flap should approximate in a manner that establishes a tension-free closure of the surgical site.
- If using ReOss® at the time of dental implant placement, additional healing time is recommended before loading dental implant.
- In the unlikely event of an adverse reaction, report immediately according to established guidelines.


⚠️ Warning

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist or physician.

Intra-Lock® International Inc.
6560 West Rogers Circle, Bldg. 24
Boca Raton FL 33487 – USA

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