

Journal section: *Implantology*Publication Types: *Review*

doi:10.4317/medoral.19306

<http://dx.doi.org/doi:10.4317/medoral.19306>

Narrow-diameter implants: Are they a predictable treatment option? A literature review

José-Luis Sierra-Sánchez ¹, Amparo Martínez-González ¹, Fernando García-Sala Bonmatí ¹, José-Félix Mañes-Ferrer ¹, Alejandro Brotons-Oliver ²

¹ Master in Advanced Oral Implantology Professor, Universidad Europea de Valencia

² Master in Advanced Oral Implantology Director, Universidad Europea de Valencia

Correspondence:

Clínica Universitaria Odontológica de la Universidad Europea de Valencia

C/ Alfambra, 4 – bajo

46009 – Valencia, Spain

jose Luis.sierra@uem.es

Sierra-Sánchez JL, Martínez-González A, García-Sala Bonmatí F, Mañes-Ferrer JF, Brotons-Oliver A. **Narrow-diameter implants: Are they a predictable treatment option? A literature review.** *Med Oral Patol Oral Cir Bucal.* 2014 Jan 1;19 (1):e74-81.

<http://www.medicinaoral.com/medoralfree01/v19i1/medoralv19i1p74.pdf>

Received: 07/05/2013
Accepted: 28/08/2013

Article Number: 19306 <http://www.medicinaoral.com/>
© Medicina Oral S. L. C.I.F. B 96689336 - pISSN 1698-4447 - eISSN: 1698-6946
eMail: medicina@medicinaoral.com

Indexed in:

Science Citation Index Expanded
Journal Citation Reports
Index Medicus, MEDLINE, PubMed
Scopus, Embase and Emcare
Índice Médico Español

Abstract

Objective: To evaluate the predictability of narrow-diameter implants as a treatment option in routine clinical practice. A literature review was performed of studies reporting clinical results obtained with these implants. Survival rates, peri-implant bone loss and related complications were evaluated. The working hypothesis was that narrow-diameter implants offer clinical results similar to those obtained with implants of greater diameter.

Material and Methods: A Medline-PubMed search covering the period between 2002 and 2012 was carried out. Studies published in English and with a follow-up period of at least 12 months were considered for inclusion. A manual search was also conducted in different journals with an important impact factor.

Results: Twenty-one studies meeting the screening criteria were included in the literature review. A total of 2980 narrow-diameter implants placed in 1607 patients were analyzed.

Conclusions: The results obtained from the literature indicate that narrow-diameter implants are a predictable treatment option, since they afford clinical results comparable to those obtained with implants of greater diameter.

Key words: *Narrow implant, survival rate, peri-implant bone loss, related complications.*

Introduction

Treatment with dental implants offers a predictable solution for most situations seen in routine clinical practice. However, bone availability is often a limiting factor in planning our treatments.

A number of surgical techniques allow us to increase the available bone width, such as for example expan-

sion with osteotomes (1), guided bone regeneration (2), autologous bone grafts (3), crestal expansion techniques (4), and osteogenic distraction (5,6).

Although these procedures offer good results in implantology, they are not without complications (7,8). The associated inconveniences are increased morbidity, longer healing times, and infection secondary to wound dehisc-

cence or membrane exposure (9,10). In patients with deficient crest width, the utilization of narrow-diameter implants therefore constitutes a technically more simple treatment alternative.

The definition of a narrow-diameter implant is subject to controversy. Although no universally accepted classification of implant diameters has been established to date, a narrow-diameter implant is generally taken to have a diameter of ≥ 3.0 mm and ≤ 3.5 mm. Some recent studies with narrow-diameter implants have reported implant success and survival rates similar to those obtained with greater diameter implants (11,12). However, few studies with prolonged periods of follow-up evaluating the predictability of these implants have been published.

The objective of the present literature review is to evaluate the predictability of narrow-diameter implants as an alternative to other technically more complex procedures, based on the survival rates, changes in peri-implant bone height and related complications.

Material and Methods

- Search strategy

A Medline-PubMed search was conducted of studies published in English between January 2002 and June 2012 (both included), using the following MESH terms: "narrow implant", "survival rate", "peri-implant bone loss", "related complications". In order to minimize search bias, we also conducted a manual search of relevant articles published in four major implant journals with an important impact factor (International Journal of Oral & Maxillofacial Implants, Clinical Oral Implants Research, Journal of Periodontology, and Clinical Implant Dentistry and Related Research). The electronic and manual searches yielded fifty-one potentially relevant studies, based on the review of the corresponding abstracts. Following the full-text evaluation of these publications, only twenty-one studies were found to meet the inclusion criteria and were finally included in the review.

- Selection of studies and inclusion criteria

A single reviewer carried out the search. The variables of interest were implant survival, changes in peri-implant bone height, and related complications. Implant survival was defined as implant persistence in the mouth at the time of evaluation.

The studies included in the review were required to meet the following criteria:

- o Full-text articles published in English in indexed journals between January 2002 and June 2012 (both included).
- o Presentation of clinical results with implants of diameter ≤ 3.5 mm and \geq of 3.0 mm.
- o Systematic reviews, randomized clinical trials, and prospective or retrospective human cohort studies.
- o A duration of follow-up of at least 12 months.

- Data extraction

All of the included studies were independently reviewed and analyzed. We collected data referred to the design of the studies (type of study, duration of follow-up, number of patients, number of implants, and type of edentulism) and the characteristics of treatment (type of implant, surgical technique employed, location of the implants, and type of prosthetic restoration). We also analyzed the variables reflecting the results of treatment (survival rate, peri-implant bone loss, associated complications [biological, prosthetic and aesthetic] and the results obtained with immediate loading).

Results

- Variables associated to the study design (Table 1)

The review included four randomized clinical trials (13-16), ten prospective studies (17-26) and seven retrospective studies (12,27-32). The follow-up periods ranged from 12 months to 12 years.

The twenty-one reviewed studies included a total of 1607 patients. Five studies involved over 100 patients (12,19,27,28,30). The patient age ranged from 13 to 87 years. Medically compromised subjects were excluded in all the studies. A total of 2980 implants were included in the review.

Regarding the type of edentulism, three studies included only totally edentulous individuals (16,21,24), five studies included both partially and totally edentulous patients (13,19,22,30,31), seven studies included partially edentulous cases (12,18,20,26-29), and six studies included only patients with a single missing tooth (13-15,23,25,32).

- Variables associated to the characteristics of treatment (Table 2)

Use was made of implants with different designs and surface treatments. The diameter of the implants ranged between 3.0 and 3.5 mm, and the length between 8 and 18 mm. Apart from the study published by Comfort et al. (17), in which machined screw-shaped narrow implants (Brånemark System) were used, twenty studies described results corresponding to implants with different surface treatments (12-16,18-32). Two of these twenty studies used both machined implants and implants with surface treatments (27,28).

Regarding the location of the implants, nine studies presented results with implants placed only in the anterior zone (13-16,21,23,25,29,32), two studies reported results with implants positioned only in the posterior zone (20,27), and ten studies presented results corresponding to both the anterior and posterior zones (12,17-19,22,24,26,28,30,31).

Three studies presented results corresponding to implants positioned only in the mandible (16,21,32), another three studies included implants only in the upper maxilla (14,24,25), and fifteen studies included implants

Table 1. Variables associated to the study design () Mean follow-up, TE: totally edentulous, PE: partially edentulous, P-TE: partially & totally edentulous, ST: single-tooth.

AUTHOR YEAR	TYPE OF STUDY	FOLLOW-UP PERIOD	NUMBER PATIENTS / IMPLANTS	TYPE OF EDENTULISM
EL SHEIKH ET AL. 2012 (26)	PROSPECTIVE	24 MONTHS	20 / 30	TE
BARTER ET AL. 2012 (14)	PROSPECTIVE	24 MONTHS	22 / 22	PE
ZEMBIC ET AL. 2012 (20)	RANDOMIZED CLINICAL TRIAL	12 MONTHS	44 / 54	ST
LEE ET AL. 2012 (27)	RETROSPECTIVE	1-12 (4.9) YEARS	338 / 541	PE
CHIAPASCO ET AL. 2012 (28)	PROSPECTIVE	24 MONTHS	18 / 51	P-TE
AL-NAWAS ET AL. 2012 (21)	RANDOMIZED CLINICAL TRIAL	12 MONTHS	89/178	TE
GALINDO-MORENO ET AL. 2011 (15)	RANDOMIZED CLINICAL TRIAL	12 MONTHS	69 / 97	ST
MALO Y DE ARAUJO NOBRE 2011 (16)	RETROSPECTIVE	1-11 (5) YEARS	147 / 247	PE
SOHN ET AL. 2011 (29)	RETROSPECTIVE	23 ± 4.3 MONTHS	32 / 62	PE
ARISAN ET AL. 2010 (30)	RETROSPECTIVE	9.1 YEARS	139 / 316	P-TE
DEGIDI ET AL 2009 (19)	RANDOMIZED CLINICAL TRIAL	3 YEARS	60 / 60	ST
DEGIDI ET AL. 2009 (18)	PROSPECTIVE	4 YEARS	40 / 93	PE
FRANCO ET AL. 2009 (31)	RETROSPECTIVE	25 MONTHS	36 / 94	P-TE
DEGIDI ET AL. 2008 (12)	RETROSPECTIVE	20 MONTHS	237 / 510	PE
REDDY ET AL. 2008 (32)	PROSPECTIVE	12 MONTHS	12 / 31	ST
VELTRI ET AL. 2008 (33)	PROSPECTIVE	12 MONTHS	12 / 73	TE
ZARONE ET AL. 2006 (34)	PROSPECTIVE	24-39 MONTHS	30 / 34	ST
ROMEO ET AL. 2006 (35)	PROSPECTIVE	7 YEARS	68 / 122	PE
CORDARO ET AL. 2006 (36)	RETROSPECTIVE	23 MONTHS	31 / 44	ST
COMFORT ET AL.2005 (13)	PROSPECTIVE	5 YEARS	9 / 23	P-TE
ZINSLI ET AL. 2004 (17)	PROSPECTIVE	6 YEARS	154 / 298	P-TE

in both the upper maxilla and mandible (12-14,17-20,22,23,26-31).

Regarding the surgical technique employed, a full-thickness flap was raised for implant placement in all the studies (12-32). Nine studies used a single-step surgical protocol (12-15,20-22,26,27), five studies used a two-step protocol (17,19,24,25,31), two studies included both single and two-step surgical protocols (18,30), and three studies performed surgery involving the raising

of a full-thickness flap without offering further details (16,28,29). Only three studies presented results comparing implants placed with a flapless technique versus the raising of a full-thickness flap (15,23,29).

Six of the reviewed studies included only single-tooth cases (13,14,15,23,25,26). Regarding the rest of the studies, three of them included only fixed partial prostheses (12,18,20), two studies presented only results corresponding to mandibular overdentures (16,21), and one

Table 2. Variables associated to the characteristics of the treatment. OD: overdenture FPP: fixed partial prostheses ANT & POST: anterior & posterior FTF: full-thickness flap.

AUTHOR YEAR	TYPE OF IMPLANT	LOCATION	SURGICAL TECHNIQUE	TYPE OF RESTORATION
EL SHEIKH ET AL. 2012 (26)	∅ 3,3 10, 12, 14 mm	MANDIBLE ANTERIOR	FTF SINGLE-STEP	OD
BARTER ET AL. 2012 (14)	∅ 3,3 10, 12 mm	BOTH JAWS ANT & POST	FTF SINGLE-STEP TWO-STEP	FPP
ZEMBIC ET AL. 2012 (20)	∅ 3,0 13, 15 mm	BOTH JAWS ANTERIOR	FTF FLAPLESS	SINGLE-TOOTH
LEE ET AL. 2012 (27)	∅3,3 ∅3,4 x3,5 8, 9, 10, 11, 11,5, 12, 13, 14, 15 mm	BOTH JAWS ANT & POST	FTF	SINGLE-TOOTH FPP
CHIAPASCO ET AL. 2012 (28)	∅ 3,3 8, 10, 12 mm	BOTH JAWS ANT & POST	FTF SINGLE-STEP	SINGLE-TOOTH FPP OD
AL-NAWAS ET AL. 2012 (21)	∅3.3 8,9,10,12,14 mm	MANDIBLE ANTERIOR	FTF	OD
GALINDO-MORENO ET AL. 2011 (15)	∅ 3,3 11, 13, 15 mm	BOTH JAWS ANTERIOR	FTF SINGLE-STEP	SINGLE-TOOTH
MALO Y DE ARAUJO NOBRE 2011 (16)	∅ 3,3 10, 11,5, 13, 15 mm	BOTH JAWS POSTERIOR	FTF	SINGLE-TOOTH FPP
SOHN ET AL. 2011 (29)	∅ 3,0 12, 15 mm	BOTH JAWS ANTERIOR	FTF FLAPLESS	SINGLE-TOOTH FPP
ARISAN ET AL. 2010 (30)	∅ 3,3 ∅3,4 8, 9,5, 11, 13, 15	BOTH JAWS ANT & POST	FTF SINGLE-STEP TWO-STEP	FPP OD
DEGIDI ET AL 2009 (18)	∅ 3,0 13,15 mm	UPPER MAXILLA ANTERIOR	FTF	SINGLE-TOOTH

Table 2. Continue.

DEGIDI ET AL. 2009 (19)	⊙ 3,0 11, 13, 15 mm	BOTH JAWS POSTERIOR	FTF	FPP
FRANCO ET AL. 2009 (31)	⊙ 3,0- ⊙ 3,5 10,13,16 mm	BOTH JAWS ANT & POST	FTF TWO-STEP	FPP OD
DEGIDI ET AL. 2008 (12)	⊙ 3,3- ⊙ 3,5 8-18 mm	BOTH JAWS ANT & POST	FTF	FPP
REDDY ET AL. 2008 (32)	⊙ 3,0	BOTH JAWS ANTERIOR	FTF FLAPLESS	SINGLE-TOOTH
VELTRI ET AL. 2008 (33)	⊙ 3,5 9, 13, 15, 17 mm	UPPER MAXILLA ANT & POST	FTF TWO-STEP	FULL ARCH FIXED REHABILITATION
ZARONE ET AL. 2006 (34)	⊙ 3,3 10, 12, 14 mm	UPPER MAXILLA ANTERIOR	FTF TWO-STEP	SINGLE-TOOTH
ROMEO ET AL. 2006 (35)	⊙ 3,3 10,12 mm	BOTH JAWS ANT & POST	FTF SINGLE-STEP	SINGLE-TOOTH FPP
CORDARO ET AL. 2006 (36)	⊙ 3,5	MANDIBLE ANTERIOR	FTF	SINGLE-TOOTH
COMFORT ET AL. 2005 (13)	⊙ 3,3 10,13,15 mm	BOTH JAWS ANT & POST	FTF TWO-STEP	FULL ARCH FIXED REHABILITATION FPP
ZINSLI ET AL. 2004 (17)	⊙ 3,3 8,10,12 mm	BOTH JAWS ANT-POST	FTF TWO-STEP	SINGLE-TOOTH FPP OD

study included only cases of full arch fixed rehabilitations (24). On the other hand, six studies presented results corresponding to both single-tooth cases and fixed partial prostheses (26,27,28,29-31), two studies presented results corresponding to single-tooth cases, fixed partial prostheses and overdentures (19,22), and one study presented results corresponding to both full arch fixed rehabilitations and fixed partial prostheses (17).

- Variables associated to the results of treatment (Table 3)
The recorded implant survival rates were above 90% in all the studies. Six studies published a survival rate of 100% at the end of the follow-up period (14,20,22,24,29,32). The lowest survival rate (90.9%) corresponded to the study published by Barter et al. (18).

58 failures (implant loss) were recorded out of a total of 2980 implants. Only five studies (12,27,28,30,31) of-

ferred information on the length of the failed implants (43 implants). A larger number of failures were recorded with implants measuring ≤ 13 mm in length, compared with implants > 13 mm in length. The authors did not examine whether these data were statistically significant or not.

Nineteen studies measured changes in peri-implant bone height after implant loading (12-18,20-30,32). The values ranged between 0.065 mm the first year according to Galindo-Moreno et al. (13) and 1.74 mm after a follow-up period of 10 years in the study published by Malo and de Araújo (27).

In relation to the recorded complications, eleven studies (13,15-17,19,21,23,26,28,30,31) registered biological complications. The latter were mainly related to a lack of implant osseointegration and infectious problems.

Table 3. Variables associated to the results of the treatment * 1 year ** 2 years *** 3 years ▪ 5 years # 10 years N.P = NOT PUBLISHED.

AUTHOR YEAR	NUMBER OF IMPLANTS	SURVIVAL RATE	PERI-IMPLANT BONE LOSS	COMPLICATIONS
EL SHEIKH ET AL. 2012 (26)	30	98 %	0.8 **	BIOLOGICAL
BARTER ET AL. 2012 (14)	22	90.9 %	0.33 ± 0.54 **	N.P
ZEMBIC ET AL. 2012 (20)	54	98 %	1.6 mm (0.8-4.6) *	BIOLOGICAL PROSTHETIC AESTHETIC
LEE ET AL. 2012 (27)	541	98.1 %	0.07 ± 0.20 *	BIOLOGICAL PROSTHETIC
CHIAPASCO ET AL. 2012 (28)	51	100 %	0-1 mm	N.P
AL-NAWAS ET AL. 2012 (21)	178	98 %	0.3 ± 0.6 mm *	BIOLOGICAL
GALINDO-MORENO ET AL. 2011 (15)	97	95.9 %	0.065 *	BIOLOGICAL PROSTHETIC
MALO Y DE ARAUJO NOBRE 2011 (16)	247	95.1 %	1.16 mm * 1.53 mm ▪ 1.74 mm #	N.P
SOHN ET AL. 2011 (29)	62	100 %	0.53 ± 0.37 mm *	N.P
ARISAN ET AL. 2010 (30)	316	92.3 %	UPPER MAXILLA 1.32 ± 0.13 mm MANDIBLE 1.28 ± 0.3 mm	BIOLOGICAL PROSTHETIC AESTHETIC
DEGIDI ET AL 2009 (19)	60	100 %	IMMEDIATE LOADING 0.85 ± 0.71 mm DELAYED LOADING 0.75 ± 0.63 mm	PROSTHETIC
DEGIDI ET AL. 2009 (18)	93	100 %	1.16 ± 0.9 mm	PROSTHETIC
FRANCO ET AL. 2009 (31)	94	95.7 %	N.P	BIOLOGICAL
DEGIDI ET AL. 2008 (12)	510	99.4 %	IMMEDIATE LOADING 1mm DELAYED LOADING 0.9mm	N.P
REDDY ET AL. 2008 (32)	31	96.7 %	0.7 mm *	BIOLOGICAL AESTHETIC
VELTRI ET AL. 2008 (33)	73	100 %	0.3 ± 0.13 *	N.P
ZARONE ET AL. 2006 (34)	34	97,06 %	1.2 ± 0.61 mm **	AESTHETIC
ROMEO ET AL. 2006 (35)	122	98.1 %	1.5 ± 1.5	BIOLOGICAL PROSTHETIC
CORDARO ET AL. 2006 (36)	44	100 %	0.8 **	AESTHETIC
COMFORT ET AL. 2005 (13)	23	95,6 %	0.41± 0.17 *	BIOLOGICAL
ZINSLI ET AL. 2004 (17)	298	96.6 %	N.P	BIOLOGICAL PROSTHETIC IMPLANT FRACTURE

Prosthetic complications were registered in eight studies (13-15,19,20,26,28,30). Only the study published by Zinsli et al. (19) reported implant fractures. Specifically, the authors recorded two fractures after an observation period of less than two years in one case and more than six years in the other.

Other prosthetic complications of lesser importance associated to the use of narrow-diameter implants were also documented, such as screw loosening, prostheses decementation, screw fracture or prostheses fracture.

Five studies presented results referred to aesthetic complications (15,23,25,30,32). These problems mainly involved a poor aesthetic outcome of the definitive restoration or insufficient filling of the interproximal papilla.

Lastly, as regards of the results obtained with immediate loading, only five studies (12,14,15,20,27) offered information comparing an immediate loading protocol and a delayed loading protocol. In 2008, Malo and de Araujo (27) published a retrospective study on 3.3 mm implants placed in posterior areas and rehabilitated following an immediate loading protocol. The recorded implant survival rate was 95.5% after 9 years of follow-up. That same year, Degidi et al. (12) published a study comparing delayed and immediate loading in narrow-diameter implants placed in both anterior and posterior zones. The reported survival rate was 99.4% after a mean follow-up of 20 months. In 2009, the same authors (20) reported a 100% survival rate with 3.0 mm implants. In this study the authors treated patients with partially edentulous posterior areas rehabilitated with an immediate loading protocol followed-up on for 48 months. Likewise in 2009, Degidi et al. (14) published another study comparing peri-implant bone loss and probing depth between the two protocols. In this study, involving a follow-up period of 36 months, the differences were not statistically significant. In 2011, Zembic et al. (15) reported a survival rate of 98% for single 3.0 mm implants subjected to immediate loading and followed-up on for one year.

Discussion

Some studies have found the survival of narrow-diameter implants (3.0 - 3.5 mm range) to be comparable to those obtained with standard-diameter implants (12,19). This review only included four randomized clinical trials supporting such performance (13-16).

Thirteen of the twenty-one studies included in the review presented a mean follow-up period of at least 24 months (14,17-22,25-28,30,31). The patient sample was quite large and included both totally and partially edentulous individuals. Due to the great variety of implants analyzed and the high survival rates recorded, it is difficult to establish a relationship between the characteristics of the different surfaces and implant survival. Regarding length, the failure rate tended to increase when

using implants measuring ≤ 13 mm in length, compared with longer implants. The statistical significance of these results could not be analyzed.

In 2006, Cardaropoli et al. (33) reported a mean interproximal bone loss of 1.5 mm during the first year for standard-diameter implants. Other authors, based on finite elements analysis, found that an increase in implant diameter exerts a greater effect in terms of the reduction of stress transmission than an increase in implant length (34).

The possibility is therefore raised that smaller-diameter implants are associated to greater bone loss after functional loading. The results obtained in this review contradict this idea, since thirteen of the reviewed studies recorded a mean bone loss of ≤ 1 mm (12-14,16-18,21-24,28,29,32). The greatest peri-implant bone loss value was 1.74 mm, reported in a study published by Malo and de Araujo, with a follow-up period of 10 years (27).

Independently of the surgical technique used, good results were obtained with narrow-diameter implants placed in both the anterior and the posterior areas of both jaws.

Favorable results were also obtained with narrow-diameter implants supporting different types of prosthetic restorations, some of which implied increased biomechanical demands upon the implants.

Although some authors have related the use of narrow-diameter implants with an increased risk of implant fracture (35,36), no such association was observed in this review. Indeed, only two fractures were recorded out of a total of 2980 implants. In any case, these results are probably conditioned by the study designs involved.

Since the need for bone grafts or regeneration techniques was obviated, the incidence of biological complications was relatively low. The recorded prosthetic complications were generally few and easy to resolve. Very little information was obtained on the aesthetic complications associated to the use of narrow-diameter implants.

Only five of the twenty-one studies included in the review contributed information on the results obtained with narrow-diameter implants subjected to immediate loading (12,14,15,20,27). The follow-up periods in these studies were 12 months, 20 months (mean follow-up), 36 months, 48 months and 9 years. Good results (with implant survival rates of over 95%) were observed in all of them. However, since the publications were so few, further studies involving longer periods of follow-up are needed in order to confirm these results.

Different factors may have influenced the results obtained. Firstly, most of the studies excluded smokers, patients with bruxism and medically compromised subjects from the analysis of results. Secondly, bone quality - which conditions primary implant stability - and the experience of the clinician, may have exerted a decisive influence.

Conclusion

Despite the limitations inherent to reviews of this kind, the results obtained appear to confirm the idea that treatment with narrow-diameter implants offers clinical results in terms of implant survival, peri-implant bone loss and associated complications similar to those of treatment with implants of greater diameter. Further studies are needed, with longer follow-up periods, in order to confirm these conclusions.

References

- Summers R.B. A new concept in maxillary Implant surgery: the osteotome technique. *Compendium*1994;15:152-8.
- Hammerle CH, Jung RE, Feloutzis A. A systematic review of the survival of implants in bone sites augmented with barrier membranes (guided bone regeneration) in partially edentulous patients. *J Clin Periodontol* 2002;29 Suppl 3:226-31;discussion 232-33.
- McAllister BS, Haghghat K. Bone augmentation techniques. *J Periodontol* 2007;78:377-96.
- Scipioni A, Bruschi G.B, Calesini G. The edentulous ridge expansion technique: a five-year study. *Int J Periodont Rest* 1994;14:451-9.
- Ilizarov GA. The tension-stress effect on the genesis and growth of tissues. Part I. The influence of stability of fixation and soft-tissue preservation. *Clin Orthop* 1989;81:238-49.
- Ilizarov GA. The tension-stress effect on the genesis and growth of tissue. Part II. The influence of the rate and frequency of distraction. *Clin Orthop* 1989;85:239-63.
- Chiapasco M, Casentini P, Zaniboni M. Bone augmentation procedures in Implant Dentistry. *Int J Oral Maxillofac Impl* 2009;24:237-59.
- Grunder U, Wenz B, Schupbach P. Guided bone regeneration around single-tooth implants in the esthetic zone: a case series. *Int J Periodont Res* 2011;31:613-20.
- Esposito M, Grusovin M.G, Felice P, Karatzopoulos G, Worthington H.V, Coulthard P. The efficacy of horizontal and vertical bone augmentation procedures for dental implants: a Cochrane systematic review. *Eur J Oral Implantol* 2009;2:167-84.
- Chiapasco M, Zaniboni M. Clinical outcomes of GBR procedures to correct peri-implant dehiscences and fenestrations: a systematic review. *Clin Oral Impl Res* 2009;20:113-23.
- Sohrabi K, Mushantat A, Esfandiari S, Feine J. How successful are small-diameter implants? A literature review. *Clin Oral Impl Res* 2012;23:515-25.
- Degidi M, Piatelli A, Carinci F. Clinical outcome of narrow diameter implants: a retrospective study of 510 implants. *J Periodontol* 2008;79:49-54.
- Galindo-Moreno P, Nilsson P, King P, Becktor J, Speroni S, Schramm A et al. Clinical and radiographic evaluation of early loaded narrow diameter implants. 1-year follow-up. *Clin. Oral. Impl. Res* 2012;23:609-16.
- Degidi M, Nardi D, Piatelli A. Immediate versus one-stage restoration of small-diameter implants for a single missing maxillary lateral incisor: a 3-year randomized clinical trial. *J. Periodontol* 2009;9:1393-8.
- Zembic A, Johannesen L.H, Schou S, Malo P, Reichert T, Farella M et al. Immediately restored one-piece single-tooth implants with reduced diameter: one-year results of a multi-center study. *Clin. Oral. Impl. Res* 2012;23:49-54.
- Al-Nawaw B, Brägger U, Meijer H.J, Naert I, Persson R, Perucchi A et al. A double-blind randomized controlled trial (RCT) of titanium-13 zirconium versus titanium grade IV small-diameter bone level implants in edentulous mandibles - results from a 1-year observation period *Clin. Impl. Dent. Rel. Res* 2012;14:896-904.
- Comfort M.B, Chu F.C.S, Chai J, Wat P.Y.P, Chow T.W. A 5-year prospective study on small diameter screw-shaped oral implants. *J Oral Rehabil* 2005;32:341-5.
- Barter S, Stone P, Brägger U. A pilot study to evaluate the success and survival rate of titanium-zirconium implants in partially edentulous patients: results after 24 months of follow-up. *Clin. Oral. Impl. Res* 2012;23:873-81.
- Zinsli B, Sägeser T, Mericske E, Mericske-Stern R. Clinical evaluation of small-diameter ITI implants: a prospective study. *Int. J. Oral. Maxillofac. Impl* 2004;19:92-9.
- Degidi M, Nardi D, Piatelli A. Immediate restoration of small-diameter implants in cases of partial posterior edentulism: a 4-year case series. *J. Periodontol* 2009;6:1006-12.
- El-Sheikh A.M, Shihabuddin O.F, Ghoraba S.M.F. Two versus three narrow-diameter Implants with locator attachments supporting mandibular overdentures. A two-year prospective study. *J. Dent* 2012. Epub ahead of print.
- Chiapasco M, Casentini P, Zaniboni M, Corsi E, Anello T. Titanium-zirconium alloy narrow-diameter implants (Straumann Roxolid®) for the rehabilitation of horizontally deficient edentulous ridges: prospective study on 18 consecutive patients. *Clin. Oral. Impl. Res* 2012;23:1136-41.
- Reddy M.S, O'Neal J, Haigh S, Aponte-Wesson R, Geurs N.C. Initial clinical efficacy of 3-mm implants immediately placed into function in conditions of limited spacing. *Int. J. Oral. Maxillofac. Impl* 2008;23:281-8.
- Veltri M, Ferrari M, Balleri P. One-year outcome of narrow diameter blasted implants for rehabilitation of maxillas with knife-edge resorption. *Clin. Oral. Impl. Res* 2008;19:1069-73.
- Zarone F, Sorrentino R, Vaccaro F, Russo S. Prosthetic treatment of maxillary lateral incisor agenesis with osseointegrated implants: a 24-39 month prospective clinical study. *Clin. Oral. Impl. Res* 2006;17:94-101.
- Romeo E, Lops D, Amorfini L, Chiapasco M, Chisolmi M, Vogel G. Clinical and radiographic evaluation of small-diameter (3.3-mm) implants followed for 1-7 years: a longitudinal study. *Clin. Oral. Impl. Res* 2006;17:139-48.
- Maló P & de Araújo Nobre M. Implants (3.3 mm diameter) for the rehabilitation of edentulous posterior regions: a retrospective clinical study with up to 11 years of follow-up. *Clin. Implant. Dent. Relat. Res* 2011;13:95-103.
- Lee J.S, Kim H.M, Kim C.S, Choi S.H, Chai J.K, Jung U.W. Long-term retrospective study of narrow implants for fixed dental prostheses. *Clin. Oral. Impl. Res* 2012. *Clin. Oral. Impl. Res* 2013;24:847-52.
- Sohn D.S, Bae M.S, Heo J.U, Park J.S, Yea S.H, Romanos G.E. Retrospective multicenter analysis of immediate provisionalization using one-piece narrow-diameter (3.0-mm) implants. *Int. J. Oral. Maxillofac. Impl* 2011;26:163-8.
- Arisan V, Bölükbası N, Ersanli S, Özdemir T. Evaluation of 316 narrow diameter implants followed for 5-10 years: a clinical and radiographic retrospective study. *Clin. Oral. Impl. Res* 2010;21:296-307.
- Franco M, Viscioni A, Rigo L, Guidi R, Zollino I, Avantaggiato A et al. Clinical outcome of narrow diameter implants inserted into allgrafts. *J. Appl. Oral. Sci* 2009;17:301-6.
- Cordaro L, Torsello F, Mirisiola di Torresanto V, Rossini C. Retrospective evaluation of mandibular incisor replacement with narrow neck implants. *Clin. Oral. Impl. Res* 2006;17:730-5.
- Cardaropoli G, Lekholm U, Wennstrom J.L. Tissue alterations at Implant-supported single-tooth replacements: a 1-year prospective clinical study. *Clin. Oral. Impl. Res* 2006;17:165-71.
- Petrie C.S, Williams J.L. Comparative evaluation of Implant designs: influence of diameter, length and taper on strains in the alveolar crest. A three-dimensional finite-element analysis. *Clin. Oral. Impl. Res* 2005;16: 486-94.
- Schwarz MS. Mechanical complications of dental implants. *Clin. Oral. Impl. Res* 2000;11:156-8.
- Allum S.R, Tomlinson R.A, Joshi R. The impact of loads on standard diameter, small diameter and mini-implants: a comparative laboratory study. *Clin. Oral. Impl. Res* 2008;19:553-9.